

VISCOLAB 3000 LABORATORY VISCOMETER

FOR VISCOSITY MEASUREMENTS FROM 0.2 TO 20,000 CENTIPOISE



Cambridge Viscosity's most advanced sensor and electronics give the new ViscoLab 3000 power to measure an exceptional range of fluids and viscosities with one lab instrument.

Operating Principle

The ViscoLab 3000 contains a piston-style electromagnetic sensor and RTD that provides continuous viscosity, temperature and temperature compensated viscosity (TCV) readings.

Measurements can be made in any of thirteen 20:1 viscosity ranges and temperature can be displayed in °C or °F.

Integrated temperature control allows user defined measurements at any setting from slightly above ambient to 356°F (180°C) and from ambient to -40°C with the optional cooling jacket. The sensor is easily field calibrated. Multiple operating ranges can be precalibrated.

Data can be output to a PC via the RS-232 serial port for tracking, storage or graphical display. The electronics allow the user to define TCV equation values, data averaging interval, alarm points, density and much more.

Key Features

- Integrated temperature control
- Requires less than 2 ml of fluid
- Simple to use, easy to clean
- Minimum bench space required
- Displays centipoise or centistokes
- Temperature display in °F or °C
- RS232 Interface
- Easily field calibrated
- Removable jacket and cap insulation

Correlates with ASTM D445

Specifications

Viscosity	0.2 to 20,000 centipoise (cP)
Piston Ranges	0.2-2cP, 0.25-5cP, 0.5-10cP, 1-20cP, 2.5-50cP, 5-100cP, 10-200cP, 25-500cP, 50-1,000cP, 100-2,000cP, 250-5,000cP, 500-10,000cP, 1,000-20,000cP
Viscosity Accuracy ¹	±1.0% of full scale
Repeatability	±0.8% of Reading
Temperature Sensor	4 wire internal Platinum RTD
Wetted Materials	316L and 430 Stainless Steel
Maximum Particle Size	25-360 Microns (Range Dependent)
Temperature Control*	Slightly above ambient to 356°F (180°C)
Standard Cable Length	3 Feet (0.9 meters)
Power	100-240 VAC, 50-60Hz

1) ± 5% of full scale 1,000-20,000cP range *ViscoLab 3000 system only

Options

Multiple Ranges (One piston per range)
Additional Cable
Insulator Sleeve (For operating temperature above 80°C)
Insulator Cap (For operating temperature above 80°C)
Cooling Jacket (For re-circulation fluids)



rev. 4/12/2006



101 Station Landing
Medford, MA 02155 USA
Tel: 781-393-6500 Fax: 781-393-6515
www.cambridgeviscosity.com
Email: info@cambridgeviscosity.com